

Product Datasheet

Guinea pig Gastrin (GT) ELISA Kit (orb1665969)

Description	Guinea pig Gastrin (GT) ELISA Kit
Reactivity	Guinea pig
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 5.9 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.