

Product Datasheet

Porcine Glutathione S Transferase Pi (GSTp) ELISA Kit (orb1665974)

Description	Porcine Glutathione S Transferase Pi (GSTp) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5-4000 pg/mL. Sensitivity: 22.8 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.