

Product Datasheet

Porcine Glutathione Peroxidase 1 (GPX1) ELISA Kit (orb1666036)

Description	Porcine Glutathione Peroxidase 1 (GPX1) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: 0.68 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.