

Product Datasheet

Human G Protein Coupled Receptor 109B (GPR109B) ELISA Kit (orb1666043)

Description Human G Protein Coupled
Receptor 109B (GPR109B) ELISA
Kit

Reactivity Human

**Tested
Applications** ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL.
Sensitivity: 0.061 ng/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 100 µL. Sample Type:
tissue homogenates, cell lysates
or other biological fluids.

Expiration Date Please enquire.