

Product Datasheet

SC11A antibody (orb166606)

Description	Rabbit polyclonal antibody to SC11A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SEC11A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	20625
Uniprot ID	P67812
NCBI	NP_001258847.1 , NP_055115.1 , NP_001258848.1 , NP_001258851.1 , NP_001258849.1 , NP_001258850.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

