

Product Datasheet

Human Gonadotropin Releasing Hormone Receptor (GnRHR) ELISA Kit (orb1666076)

Description Human Gonadotropin Releasing Hormone Receptor (GnRHR) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL. Sensitivity: 0.07 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.

Expiration Date Please enquire.