

Product Datasheet

General Glucocorticoid (GC) ELISA Kit (orb1666098)

| | |
|----------------------------|--|
| Description | General Glucocorticoid (GC) ELISA Kit |
| Reactivity | All |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 0.156-10 ng/mL. Sensitivity: 0.059 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Expiration Date | Please enquire. |