

Product Datasheet

Mouse Glutaminase 2 (GLS2) ELISA Kit (orb1666102)

Description Mouse Glutaminase 2 (GLS2) ELISA Kit

Reactivity Mouse

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.625-40 ng/mL.
Sensitivity: 0.25 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates or other biological fluids.

Expiration Date Please enquire.