

## Product Datasheet

# Mouse Gamma-Aminobutyric Acid A Receptor Alpha 2 (gABRa2) ELISA Kit (orb1666267)

|                            |                                                                                                                                                                                                             |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Gamma-Aminobutyric Acid A Receptor Alpha 2 (gABRa2) ELISA Kit                                                                                                                                         |
| <b>Reactivity</b>          | Mouse                                                                                                                                                                                                       |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                       |
| <b>Note</b>                | For research use only                                                                                                                                                                                       |
| <b>Application notes</b>   | Detect Range: 0.156-10 ng/mL. Sensitivity: 0.051 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                             |