

## Product Datasheet

# Human Gamma-Aminobutyric Acid A Receptor Alpha 2 (gABRa2) ELISA Kit (orb1666268)

<b>Description</b>	Human Gamma-Aminobutyric Acid A Receptor Alpha 2 (gABRa2) ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.059 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
<b>Expiration Date</b>	Please enquire.