

Product Datasheet

Human Ferritin, Heavy Polypeptide (FTH) ELISA Kit (orb1666300)

Description	Human Ferritin, Heavy Polypeptide (FTH) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 12.9 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.