

Product Datasheet

General Free Thyroxine (fT4) ELISA Kit (orb1666301)

Description General Free Thyroxine (fT4) ELISA Kit

Reactivity All

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 1.4-90 pg/mL.
Sensitivity: 0.48 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Expiration Date Please enquire.