

Product Datasheet

Porcine Follicle Stimulating Hormone (FSH) ELISA Kit (orb1666307)

Description	Porcine Follicle Stimulating Hormone (FSH) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 6.25-400 mIU/mL. Sensitivity: 2.13 mIU/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.