

## Product Datasheet

# Porcine Fibroblast Growth Factor 7 (FGF7) ELISA Kit (orb1666378)

**Description** Porcine Fibroblast Growth Factor 7 (FGF7) ELISA Kit

**Reactivity** Porcine

**Tested Applications** ELISA

**Note** For research use only

**Application notes** Detect Range: 15.6-1000 pg/mL. Sensitivity: 5.9 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.

**Expiration Date** Please enquire.