

Product Datasheet

Porcine Fibroblast Growth Factor 2, Basic (FGF2) ELISA Kit (orb1666389)

Description	Porcine Fibroblast Growth Factor 2, Basic (FGF2) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 41.78 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.