

Product Datasheet

Human Fibroblast Growth Factor 17 (FGF17) ELISA Kit (orb1666410)

Description	Human Fibroblast Growth Factor 17 (FGF17) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 35.96 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.