

Product Datasheet

Human Fc Fragment Of IgG Binding Protein (FcγBP) ELISA Kit (orb1666443)

Description	Human Fc Fragment Of IgG Binding Protein (FcγBP) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.116 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.