

Product Datasheet

Human Fibulin 4 (FBLN4) ELISA Kit (orb1666455)

| | |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Human Fibulin 4 (FBLN4) ELISA Kit |
| Reactivity | Human |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 3.12-200 ng/mL. Sensitivity: 1.48 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids. |
| Expiration Date | Please enquire. |