

Product Datasheet

Human Exosome Component 2 (EXOSC2) ELISA Kit (orb1666543)

Description	Human Exosome Component 2 (EXOSC2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.061 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.