

Product Datasheet

General Epitestosterone (ET) ELISA Kit (orb1666551)

Description General Epitestosterone (ET) ELISA Kit

Reactivity All

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 31.2-2000 pg/mL.
Sensitivity: 8.33 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.

Expiration Date Please enquire.