

## Product Datasheet

# Mouse EGF Like Domain Protein, Multiple 7 (EGFL7) ELISA Kit (orb1666648)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Mouse EGF Like Domain Protein, Multiple 7 (EGFL7) ELISA Kit  |
| <b>Reactivity</b>          | Mouse  |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Detect Range: 31.2-2000 pg/mL. Sensitivity: 15.5 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.  |