

Product Datasheet

Human Endothelin 2 (EDN2) ELISA Kit (orb1666661)

Description	Human Endothelin 2 (EDN2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 11.6 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.