

Product Datasheet

Mouse Endothelin 1 (EDN1) ELISA Kit (orb1666664)

Description	Mouse Endothelin 1 (EDN1) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 7.81-500 pg/mL. Sensitivity: 2.45 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.