

Product Datasheet

Monkey Eosinophil Chemotactic Factor (ECF) ELISA Kit (orb1666669)

Description	Monkey Eosinophil Chemotactic Factor (ECF) ELISA Kit
Reactivity	Monkey
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 7.3 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.