

Product Datasheet

Equine Eosinophil Chemotactic Factor (ECF) ELISA Kit (orb1666674)

Description	Equine Eosinophil Chemotactic Factor (ECF) ELISA Kit
Reactivity	Equine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5-4000 pg/mL. Sensitivity: 26.1 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.