

Product Datasheet

Rat Big Dynorphin (Dyn) ELISA Kit (orb1666685)

| | |
|----------------------------|---|
| Description | Rat Big Dynorphin (Dyn) ELISA Kit |
| Reactivity | Rat |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 7.81-500 pg/mL. Sensitivity: 3.3 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids. |
| Expiration Date | Please enquire. |