

Product Datasheet

Human Dual Specificity Phosphatase 1 (DUSP1) ELISA Kit (orb1666695)

Description	Human Dual Specificity Phosphatase 1 (DUSP1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 12.9 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.