

Product Datasheet

Human Delta-Sleep Inducing Peptide (dSIP) ELISA Kit (orb1666698)

Description	Human Delta-Sleep Inducing Peptide (dSIP) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 125-8000 pg/mL. Sensitivity: 47 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.