

Product Datasheet

Rat Dopamine Receptor D1 (DRD1) ELISA Kit (orb1666711)

Description	Rat Dopamine Receptor D1 (DRD1) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.062 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.