

Product Datasheet

Human Dolichyl Phosphate-N-Acetylglucosaminophosphotransferase 1 (DPAGT1) ELISA Kit (orb1666724)

Description	Human Dolichyl Phosphate-N-Acetylglucosaminophosphotransferase 1 (DPAGT1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.059 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.