

## Product Datasheet

# General 1-Desamino 8D Arginine Vasopressin (DDAVP) ELISA Kit (orb1666816)

**Description** General 1-Desamino 8D Arginine Vasopressin (DDAVP) ELISA Kit

**Reactivity** All

**Tested Applications** ELISA

**Note** For research use only

**Application notes** Detect Range: 0.156-10 ng/mL. Sensitivity: 0.039 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50  $\mu$ L. Sample Type: serum, plasma or other biological fluids.

**Expiration Date** Please enquire.