


Product Datasheet

SPNXB antibody (orb166684)

Description	Rabbit polyclonal antibody to SPNXB
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SPANXB1
Form/Appearance	Supplied in PBS with 0.09% (W/V) sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Purity	This antibody is purified through a protein A column, followed by peptide affinity purification.
MW	11840 Da
Uniprot ID	Q9NS25
NCBI	NP_115850.1 , NP_663697.1 , NP_620588.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis in MDA-MB231 cell ...