

Product Datasheet

Porcine Diamine Oxidase (DAO) ELISA Kit (orb1666842)

Description Porcine Diamine Oxidase (DAO) ELISA Kit

Reactivity Porcine

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 1.56-100 ng/mL. Sensitivity: 0.67 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates or other biological fluids.

Expiration Date Please enquire.