

Product Datasheet

Rabbit D-Dimer (D2D) ELISA Kit (orb1666852)

Description	Rabbit D-Dimer (D2D) ELISA Kit
Reactivity	Rabbit
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.12-200 ng/mL. Sensitivity: 1.27 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: plasma.
Expiration Date	Please enquire.