

## Product Datasheet

# Mouse Cysteine Rich Protein, Angiogenic Inducer 61 (CYR61) ELISA Kit (orb1666859)

**Description** Mouse Cysteine Rich Protein, Angiogenic Inducer 61 (CYR61) ELISA Kit

**Reactivity** Mouse

**Tested Applications** ELISA

**Note** For research use only

**Application notes** Detect Range: 15.6-1000 pg/mL. Sensitivity: 5.4 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100  $\mu$ L. Sample Type: serum, plasma, tissue homogenates or other biological fluids.

**Expiration Date** Please enquire.