

Product Datasheet

Rat Chemokine C-X-C-Motif Receptor 3 (CXCR3) ELISA Kit (orb1666891)

Description	Rat Chemokine C-X-C-Motif Receptor 3 (CXCR3) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.059 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.