

Product Datasheet

Rat Cross Linked C-Telopeptide Of Type II Collagen (CTXII) ELISA Kit (orb1666919)

Description Rat Cross Linked C-Telopeptide Of Type II Collagen (CTXII) ELISA Kit

Reactivity Rat

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 125-8000 pg/mL. Sensitivity: 33.5 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Expiration Date Please enquire.