

Product Datasheet

Rabbit Connective Tissue Growth Factor (CTGF) ELISA Kit (orb1666946)

| | |
|----------------------------|--|
| Description | Rabbit Connective Tissue Growth Factor (CTGF) ELISA Kit |
| Reactivity | Rabbit |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 46.8-3000 pg/mL. Sensitivity: 18.54 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma or other biological fluids. |
| Expiration Date | Please enquire. |