

Product Datasheet

Porcine Cystatin 3 (CST3) ELISA Kit (orb1666964)

Description	Porcine Cystatin 3 (CST3) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.056 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, urine or other biological fluids.
Expiration Date	Please enquire.