

Product Datasheet

Human Complexin 2 (CPLX2) ELISA Kit (orb1667052)

Description	Human Complexin 2 (CPLX2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.084 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 μ L. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.