

Product Datasheet

Bovine Cluster Of Differentiation 40 Ligand (CD40L) ELISA Kit (orb1667331)

Description Bovine Cluster Of Differentiation 40 Ligand (CD40L) ELISA Kit

Reactivity Bovine

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.625-40 ng/mL. Sensitivity: 0.18 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates, cell lysates or other biological fluids.

Expiration Date Please enquire.