

Product Datasheet

Mouse Cluster Of Differentiation (CD163) ELISA Kit (orb1667351)

| | |
|----------------------------|--|
| Description | Mouse Cluster Of Differentiation (CD163) ELISA Kit |
| Reactivity | Mouse |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 78.1-5000 ng/mL. Sensitivity: 25 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids. |
| Expiration Date | Please enquire. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)