

Product Datasheet

Mouse Chemokine C-C-Motif Receptor 1 (CCR1) ELISA Kit (orb1667370)

Description	Mouse Chemokine C-C-Motif Receptor 1 (CCR1) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.057 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.