

Product Datasheet

Human Cholecystokinin 8, Octapeptide (CCK8) ELISA Kit (orb1667401)

Description	Human Cholecystokinin 8, Octapeptide (CCK8) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 4.46 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.