

Product Datasheet

Guinea pig Cholecystokinin 8, Octapeptide (CCK8) ELISA Kit (orb1667402)

Description	Guinea pig Cholecystokinin 8, Octapeptide (CCK8) ELISA Kit
Reactivity	Guinea pig
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 5.28 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.