

## Product Datasheet

### Human Caveolin 1 (CAV1) ELISA Kit (orb1667424)

<b>Description</b>	Human Caveolin 1 (CAV1) ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.234-15 ng/mL. Sensitivity: 0.07 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Expiration Date</b>	Please enquire.