

Product Datasheet

Porcine Bone Morphogenetic Protein 4 (BMP4) ELISA Kit (orb1667614)

Description	Porcine Bone Morphogenetic Protein 4 (BMP4) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.067 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.