

Product Datasheet

Human Aspartate Beta Hydroxylase (ASPH) ELISA Kit (orb1667761)

| | |
|----------------------------|---|
| Description | Human Aspartate Beta Hydroxylase (ASPH) ELISA Kit |
| Reactivity | Human |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 0.156-10 ng/mL. Sensitivity: 0.057 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates or other biological fluids. |
| Expiration Date | Please enquire. |