

Product Datasheet

Human Aquaporin 2, Collecting Duct (AQP2) ELISA Kit (orb1667804)

Description	Human Aquaporin 2, Collecting Duct (AQP2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.122 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: tissue homogenates, urine or other biological fluids.
Expiration Date	Please enquire.