

Product Datasheet

Porcine Apolipoprotein A2 (APOA2) ELISA Kit (orb1667845)

| | |
|----------------------------|--|
| Description | Porcine Apolipoprotein A2 (APOA2) ELISA Kit |
| Reactivity | Porcine |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 4.68-300 µg/mL. Sensitivity: 1.2 µg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids. |
| Expiration Date | Please enquire. |